


BONASTENT® Tracheal/Bronchial

Coding: B-BONASTENT / T-Tracheal / B-Bronchial



Order Number	Stent		Delivery Device	
	Diameter [mm]	Total Length [mm]	Diameter [mm]	Usable Length [mm]
BTB-100209	10	20	2.66 [8Fr]	900
BTB-100309	10	30	2.66 [8Fr]	900
BTB-100409	10	40	2.66 [8Fr]	900
BTB-100509	10	50	2.66 [8Fr]	900
BTB-100609	10	60	2.66 [8Fr]	900
BTB-100709	10	70	2.66 [8Fr]	900
BTB-100809	10	80	2.66 [8Fr]	900
BTB-120209	12	20	3.00 [9Fr]	900
BTB-120309	12	30	3.00 [9Fr]	900
BTB-120409	12	40	3.00 [9Fr]	900
BTB-120509	12	50	3.00 [9Fr]	900
BTB-120609	12	60	3.00 [9Fr]	900
BTB-120709	12	70	3.00 [9Fr]	900
BTB-120809	12	80	3.00 [9Fr]	900
BTB-140209	14	20	3.35 [10Fr]	900
BTB-140309	14	30	3.35 [10Fr]	900
BTB-140409	14	40	3.35 [10Fr]	900
BTB-140509	14	50	3.35 [10Fr]	900
BTB-140609	14	60	3.35 [10Fr]	900
BTB-140709	14	70	3.35 [10Fr]	900
BTB-140809	14	80	3.35 [10Fr]	900
BTB-160209	16	20	4.00 [12Fr]	900
BTB-160309	16	30	4.00 [12Fr]	900
BTB-160409	16	40	4.00 [12Fr]	900
BTB-160509	16	50	4.00 [12Fr]	900
BTB-160609	16	60	4.00 [12Fr]	900
BTB-160709	16	70	4.00 [12Fr]	900
BTB-160809	16	80	4.00 [12Fr]	900
BTB-180209	18	20	5.20 [16Fr]	900
BTB-180309	18	30	5.20 [16Fr]	900
BTB-180409	18	40	5.20 [16Fr]	900
BTB-180509	18	50	5.20 [16Fr]	900
BTB-180609	18	60	5.20 [16Fr]	900
BTB-180709	18	70	5.20 [16Fr]	900
BTB-180809	18	80	5.20 [16Fr]	900
BTB-200309	20	30	5.20 [16Fr]	900
BTB-200409	20	40	5.20 [16Fr]	900
BTB-200509	20	50	5.20 [16Fr]	900
BTB-200609	20	60	5.20 [16Fr]	900
BTB-200709	20	70	5.20 [16Fr]	900
BTB-200809	20	80	5.20 [16Fr]	900